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Bre-Gel BG2 opak

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Instructions for use

Bre-Gel BG2 opaque

Before using this product, please read these instructions for use carefully!

1. Indication range

Bre-Gel BG2 opaque is a reversible duplicating gel for the fabrication of duplicate models in CoCr casting and for the use in the casting resin technique. The constant quality of Bre-Gel BG2 opaque ensures reuseability for a large number of additional models. Excellent flow characteristics and uniform solidification as well as the high stability allow precise reproduction of details of attachments and interlocks. Even deep undercuts can be exactly reproduced thanks to the high elasticity (tensile strength). The opaque color allows perfect control of the gel impression.

2. Contraindication

The duplicate mould made from Bre-Gel BG 1 can not be filled (poured) with stone.

3. Processing

3.1.1 Melting in the water bath

Cut Bre-Gel BG2 opaque into small pieces and melt in a jar in the water bath on a heating device (electric hotplate or similar device). During this process, Bre-Gel BG2 opaque must be repeatedly stirred. During heating of Bre-Gel BG2 opaque, the temperature must not exceed 93° C. The gel must not be boiled at temperatures exceeding 95° C. Careful melting will increase the service life of Bre-Gel BG2 opaque considerably.

3.1.2 Melting in the automatic duplicating unit

Cut Bre-Gel BG2 opaque into somewhat bigger pieces and fill into the duplicating unit. The heating temperature of the unit must be adjusted to 90° C. The ideal working temperature is between 40° and 45° C. Please adhere to the operating instructions for the respective duplicating unit. After heating and subsequent cooling down to the working temperature the gel is ready to use.

3.1.3 Melting in the microwave unit

Cut Bre-Gel BG2 opaque into small pieces and fill into a glass or plastic cup. Place the cup into the microwave unit and melt for 2 minutes at 900 watts. Take out the glass or plastic cup, stir Bre-Gel BG2 opaque and melt again for 1 minute a 600 watts. Immediately after Bre-Gel BG2 opaque starts to bubble, stop the heating process and remove Bre-Gel BG2 opaque from the microwave unit.

Bre-Gel BG2 opaque must not be boiled (from 95° C). Careful melting will increase the service life of Bre-Gel BG2 opaque considerably.

3.2 Duplicating

The optimum working temperature is 40° to 45° C. During cooling of Bre-Gel BG2 opaque, the material must be stirred frequently to avoid the formation of lumps. Gelation time is approx. 40 minutes at air. Then the master models must be devested and Brevest Exakta M, (REF 570 OXM2 0) is poured into the duplicating mould. Remove the investment material model from the duplicating mould.

4. Reuse and disposal

Clean the duplicating mould under running water. If Bre-Gel BG2 opaque is to be reused, it must be ensured that only clean gel material is heated again. In case of insufficient elasticity or deviations in color or changes in consistency, Bre-Gel BG2 opaque must be disposed of together with regular household garbage.

5. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore, we recommend reading the instructions for use again before using a new package. Any recommendations concerning the use are based on bredent's experience and tests and can only be considered to provide approximate values. Our products are subject to continuous further development. Therefore, we reserve the right to perform changes in the design and composition.

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